

PROTOCOL 1B

TALE/N assembly into DIRECT vectors

Time needed to complete this protocol: **6 days (2 cloning steps)**

Time needed to obtain the plant transformation vector with complete TALE/Ns: **9 days (3 cloning steps)**

Vectors compatible with this protocol: **pDIRECT_21J, pDIRECT_22J, pDIRECT_23J, pDIRECT_25K, pDIRECT_26K, pDIRECT_37J, pDIRECT_38J, pDIRECT_39J**

Summary: The first TALE will be cloned into the Esp3I sites of the pDIRECT vector, replacing the *ccdB* gene. After clones containing the first TALE are identified and DNA prepared, the second TALE is cloned into the BsaI sites of the pDIRECT-derived vector, replacing the *LacZ* gene. Correct clones are ready for plant transformation after sequencing. See also PROTOCOL 1 description.

Enzymes:

- **Esp3I**
- **BsaI**
- **T4 DNA ligase + buffer**
- **DNA polymerase (for colony PCR)**

1. Start with the repeat arrays for TALE/N 1 cloned into pFUS_A and pFUS_B plasmids (Day 3 in Cermak et al. (2015). TALE/N 2 will be assembled in another cloning step.
2. Prepare the following reaction:
 - a) 150 ng of each array vector for TALE/N 1 (pFUS_A with ten repeats cloned and pFUS_B with ten or fewer)
 - b) 150 ng of the last repeat vector
 - c) 75 ng of the pDIRECT vector
 - d) 1 µl Esp3I
 - e) 1 µl T4 DNA ligase
 - f) 2 µl 10X T4 DNA ligase buffer
 - g) H₂O up to 20 µl total reaction volume
3. Place the Golden Gate reactions in a PCR machine and run the following cycle: 10x (37°C/5min + 16°C/10min) + 37°C/15min + 80°C/5min.
4. Transform 5 µl of each Golden Gate reaction into *E. coli* (DH5α or similar, but sensitive to the presence of the *ccdB* gene) and plate on LB + 50mg/L kanamycin. Note, if X-gal is included in this step, all colonies (correct and incorrect) should be blue.

5. Correct clones can be identified via colony PCR, using primers TAL2Aseq_F1 and TAL_R2 (see the table below for primer sequences)
6. Isolate plasmid DNA for one correct clone.
7. Prepare the following reaction:
 - d) 150 ng of each array vector for TALE/N 2 (pFUS_A with ten repeats cloned and pFUS_B with ten or fewer)
 - e) 150 ng of the last repeat vector
 - f) 75 ng of the vector from step 6
 - d) 0.5 µl BsaI
 - e) 0.5 µl Esp3I
 - f) 1 µl T4 DNA ligase
 - g) 2 µl 10X T4 DNA ligase buffer
 - h) H₂O up to 20 µl total reaction volume
8. Place the Golden Gate reactions in a PCR machine and run the following cycle: 10x (37°C/5min + 16°C/10min) + 37°C/15min + 80°C/5min.
9. Transform 5 µl of each Golden Gate reaction into *E. coli* (DH5α or similar) and plate on LB + 50mg/L kanamycin + 32mg/L X-gal.
10. Correct clones can be identified via PCR on white colonies, using primers TALSeq_5-1 and TAL2A_SeqR2 (see the table below for primer sequences)
11. Isolate the plasmid DNA for one correct clone (can be sequenced using primers TAL2Aseq_F1 and TAL_R2 for TALEN 1 and TALseq_5-1 and TAL2Aseq_R2 for TALEN 2).

Primers for colony PCR and sequencing (5' to 3')

TALE/N 1	TAL2Aseq_F1	TGCCAGCAGGAATGCGTTA
	TAL_R2	GGCGACGAGGTGGTCGTTGG
TALE/N 2	TALseq_5-1	CATCGCGCAATGCACTGAC
	TAL2Aseq_R2	CGCGAGTGCCACCAGATGA

REFERENCES

Cermak, T., Starker, C.G., and Voytas, D.F. (2015). Efficient design and assembly of custom TALENS using the golden gate platform. *Methods Mol. Biol.* **1239**: 133–159.