



OsU3 promoter: TGGCXXXXXXXXXXXXXXXXXXXXX  
 |||||  
 YYYYYYYYYYYYYYYYYYYYCAA

OsU6 promoter: TTGTXXXXXXXXXXXXXXXXXXXXX  
 |||||  
 YYYYYYYYYYYYYYYYYYYYCAA

2. Phosphorylate the oligonucleotides. (Alternatively, phosphorylated oligos may be purchased. Phosphorylating them yourself is efficient enough and much more cost effective.) Prepare the following reaction:
  - a) 3  $\mu$ l 100  $\mu$ M sense gRNA oligonucleotide
  - b) 3  $\mu$ l 100  $\mu$ M antisense gRNA oligonucleotide
  - c) 3  $\mu$ l T4 DNA ligase buffer (contains ATP)
  - d) 2  $\mu$ l T4 polynucleotide kinase
  - e) 19  $\mu$ l H<sub>2</sub>O
3. Incubate 1 hour at 37°C.
4. Denature and gradually cool down the DNA. This can be done in a PCR machine using the following program: 95°C/5 min + ramping down to 85°C at -2°C/second + ramping down to 25°C at -0.1°C/second + 4°C hold **OR** by boiling the reaction in a water bath for 2 minutes and letting it cool down gradually.
5. Dilute the reaction 25 times (1  $\mu$ l oligo mixture + 24  $\mu$ l H<sub>2</sub>O)
6. Setup a Golden Gate reaction:
  - a) 50 ng of selected module B or C plasmid
  - b) 1  $\mu$ l 25x diluted annealed oligonucleotides
  - c) 0.5  $\mu$ l Esp3I (select modules use BsaI instead, check the plasmid description and/or map)
  - d) 2  $\mu$ l 10X T4 DNA ligase buffer
  - e) 1  $\mu$ l T4 DNA ligase
  - f) H<sub>2</sub>O up to 20  $\mu$ l
7. Place the Golden Gate reaction in a PCR machine and run the following cycle: 37°C/5min + 16°C/10min + 37°C/15min + 80°C/5min.
8. Transform 5  $\mu$ l of the Golden Gate reaction into *E. coli* (DH5 $\alpha$  or similar, but sensitive to the presence of ccdB gene) and plate on LB + 50mg/L ampicillin/carbenicillin.
9. Correct clones can be identified via colony PCR using the sense gRNA oligonucleotide as the forward primer and primer ZY015F (see the table below for primer sequence) as the reverse primer. However, this is usually not necessary thanks to the high cloning efficiency.

10. Isolate the plasmid DNA for one correct clone (can be sequenced using the ZY015F primer). Note that module C plasmids can be reused in another cloning step to insert a gene targeting donor following **PROTOCOL 4**.
11. Assemble gRNAs in modules B and/or C along with a Cas9 cassette in module A into selected transformation backbone using **PROTOCOL 5**.

#### Oligonucleotides in 5' to 3' orientation

AtU6 sense gRNA oligo	GATTXXXXXXXXXXXXXXXXXXXXX
At7SL sense gRNA oligo	GTACXXXXXXXXXXXXXXXXXXXXX
TaU3 sense gRNA oligo	AAGCXXXXXXXXXXXXXXXXXXXXX
TaU6 sense gRNA oligo	ACTTXXXXXXXXXXXXXXXXXXXXX
OsU3 sense gRNA oligo	TGGCXXXXXXXXXXXXXXXXXXXXX
OsU6 sense gRNA oligo	TTGTXXXXXXXXXXXXXXXXXXXXX
antisense gRNA oligo (all)	AAACYYYYYYYYYYYYYYYYYYY
ZY015F	GGAATAAGGGCGACACGGAAATG