

## PROTOCOL 4

### Assembly of DNA donor templates for gene targeting into module C vectors

Time needed to complete this protocol: **3 days (1 cloning step)**

Time needed to obtain the plant transformation vector with a donor template and SSN expression cassette: **6 days (2 cloning steps) or 9 days (3 cloning steps) if gRNA spacers are cloned in module C**

Vectors compatible with this protocol: **all module C plasmids**

**Summary:** The donor template for gene targeting will be assembled from one or more PCR fragments into the Bael site of the pMOD\_C vector by Gibson assembly. Correct clones are ready for assembly into transformation backbones (PROTOCOL 5). See also PROTOCOL 4 description.

Enzymes:

- proof reading DNA polymerase
- Bael
- 1.33x Gibson assembly master mix (make as 15µl aliquots, store at -20°C):
  - 2 µl T5 exonuclease (1 U/µl)
  - 6.25 µl Phusion DNA polymerase (2 U/µl)
  - 50 µl Taq DNA ligase (40 U/µl)
  - 100 µl 5x isothermal buffer (25% PEG-8000, 500 mM Tris-HCl pH7.5, 50 mM MgCl<sub>2</sub>, 50 mM DTT, 1 mM dATP, 1 mM dCTP, 1 mM dGTP, 1 mM dTTP, 5 mM NAD)
  - H<sub>2</sub>O to 375 µl

1. Design primers as follows:

oBael\_FWD (sequence in blue overlaps with Bael digested module C vector)

CGCGTAGTCCTCGGTAxxxxxxxxxxxxxxxxxxxx – replace red Xs with a primer binding to the 5' end of the left homology arm

Reverse primer (REV1) to be used with oBael\_FWD that overlaps with the forward primer (FWD2) for the next fragment (genomic or insertion sequence). First 8 bp of the FWD2 primer is added to the 5' end of the REV1 primer in reverse complement and vice versa to create a 16 bp overlap. T<sub>m</sub> of the overlap should be at least 48°C (add 2°C for each A-T pair and 4°C for each G-C pair to calculate the T<sub>m</sub>). If T<sub>m</sub> is below 48°C, the overlap can be extended until the T<sub>m</sub> is at least 48°C. Sequence modifications (base substitutions for allele replacement or AarI removal) can be in the overlapping region of both primers. Additional primers are designed in the same way for each other fragment (e.g. REV2 pairing with FWD2 and overlapping with FWD3, REV3 pairing with FWD3 and overlapping with FWD4, etc.). Last FWD primer will pair with oBael\_REV.



12. Assemble module C plasmid containing the donor template along with selected modules A and B into selected transformation backbone using **PROTOCOL 5**.

**Primer for colony PCR and sequencing (5' to 3')**

ZY015F	GGAATAAGGGCGACACGGAAATG
--------	-------------------------